

Anti-ADCY5 antibody (STJ502057)

STJ502057

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADCY5 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:200-1:1, 000 ICC: 1:200-1:1, 000 IF: 1:200-1:1, 000 CFM: 1:200-1:1, 000
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64532
Gene Symbol	Adcy5
Uniprot ID	ADCY5_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on rat Particulate Adenylyl Cyclase 5 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunohistochemistry with Anti-ADCY5 Antibody (STJ502057) on rat striatal neurons.