

Anti-ADCY9 antibody (STJ502061)

STJ502061

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADCY9 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-1.20 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115
Gene Symbol	ADCY9
Uniprot ID	ADCY9_HUMAN
Immunogen	Adenylate cyclase type 9 recombinant protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-ADCY9 Antibody (STJ502061) with AC9 peptide on at 1:500 dilution in buffer.