

Anti-ADCYAP1 antibody (125-176) {FITC} (STJ502082)

STJ502082

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADCYAP1 (125-176) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116
Gene Symbol	ADCYAP1
Uniprot ID	PACA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 125-176 on human Pituitary adenylate cyclase-activating polypeptide protein.
Immunogen Region	125-176
Specificity	
Immunogen Sequence	