

Anti-Phospho-alpha Actinin 4-S159 antibody (140-190) {Biotin} (STJ502087)

STJ502087

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-alpha Actinin 4-S159 (140-190) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.62 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81
Gene Symbol	ACTN4
Uniprot ID	ACTN4_HUMAN
Immunogen	Synthetic phospho-peptide, corresponding to a region within amino acid regions 140-190 of Human alpha Actinin 4 conjugated to a carrier protein.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	