

Anti-PAX1 antibody (210-260) {FITC} (STJ502124)

STJ502124

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PAX1 (210-260) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-0.95 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5075
Gene Symbol	PAX1
Uniprot ID	PAX1_HUMAN
Immunogen	Synthetic peptide corresponding to unique sequence within amino acid region 210-260 on human PAX1 protein.
Immunogen Region	210-260
Specificity	
Immunogen Sequence	