

## Anti-PAX2 isoform d + e antibody (367-417) {FITC} (STJ502145)

STJ502145

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-PAX2 isoform d + e (367-417) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.75-0.95 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5076</a>
<b>Gene Symbol</b>	<a href="#">PAX2</a>
<b>Uniprot ID</b>	<a href="#">PAX2_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 367-417 on human PAX2 isoform d protein.
<b>Immunogen Region</b>	367-417
<b>Specificity</b>	This antibody will detect Pax-2 isoform d and Pax-2 isoform e. Will not detect other Pax-2 isoforms.
<b>Immunogen Sequence</b>	