

Anti-PAX2 isoform a + e antibody (200-220) {Biotin} (STJ502147)

STJ502147

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PAX2 isoform a + e (200-220) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-0.95 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5076
Gene Symbol	PAX2
Uniprot ID	PAX2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-220 on human PAX2 isoform d protein.
Immunogen Region	200-220
Specificity	This antibody will detect Pax-2 isoform a and Pax-2 isoform e. Will not detect other Pax-2 isoforms.
Immunogen Sequence	