

Anti-PAX5 antibody (250-300) {FITC} (STJ502163)

STJ502163

GENERAL INFORMATION

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-PAX5 (250-300) is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|-----------------------------------------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.65-0.75 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---------------------------------------------------------------------------------|
| Gene ID | 5079 |
| Gene Symbol | PAX5 |
| Uniprot ID | PAX5_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 250-300 on human PAX5 protein. |
| Immunogen Region | 250-300 |
| Specificity | |
| Immunogen Sequence | |