

Anti-Phospho-CACNB2 antibody (600-655) {Biotin} (STJ502191)

STJ502191

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-CACNB2 (600-655) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116600
Gene Symbol	Cacnb2
Uniprot ID	CACB2_RAT
Immunogen	Synthetic peptide taken within amino acid region 600-655 on rat CACNB2 protein.
Immunogen Region	600-655
Specificity	
Immunogen Sequence	