

Anti-Phospho-Chimaerin beta 2 antibody {Biotin} (STJ502194)

STJ502194

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Chimaerin beta 2 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.63 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 0 IHC: 0 ICC: 0 CFM: 0
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1123
Gene Symbol	CHN1
Uniprot ID	CHIN_HUMAN
Immunogen	Phosphorylated synthetic peptide corresponding to unique amino acid sequence on beta 2 Chimaerin protein.
Immunogen Region	
Specificity	
Immunogen Sequence	