

Anti-Phospho-NT5C2 antibody {FITC} (STJ502198)

STJ502198

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Phospho-NT5C2 is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA, WB |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.75-0.85 µg/µl |
| Conjugation | FITC |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 22978 |
| Gene Symbol | NT5C2 |
| Uniprot ID | 5NTC_HUMAN |
| Immunogen | Synthetic peptide corresponding to unique amino acid sequences on cNII protein with phosphorylated Serine 7. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |