

Anti-PDE2A antibody (1-50) (STJ502259)

STJ502259

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE2A (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5138
Gene Symbol	PDE2A
Uniprot ID	PDE2A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human PDE2A4 protein.
Immunogen Region	1-50
Specificity	This antibody will detect PDE2A4 only.
Immunogen Sequence	



Western blot with Anti-PDE2A Antibody (STJ502259) and PDE2A4 peptide. Antibody at 1:500 dilution in buffer.