

Anti-PDE3B antibody (1060-1110) (STJ502269)

STJ502269

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Goat polyclonal antibody anti-PDE3B (1060-1110) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IP, WB |
| Host/Source | Goat |
| Reactivity | Chicken, Human, Monkey, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55 µg/µl |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:500 ICC: 1:500 IF: 1:500 CFM: 1:500 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 29516 |
| Gene Symbol | Pde3b |
| Uniprot ID | PDE3B_RAT |
| Immunogen | Synthetic peptide taken within amino acid region 1060-1110 on human PDE3B protein. |
| Immunogen Region | 1060-1110 |
| Specificity | This antibody does not cross react with PDE4A variants nor with any other PDE family members. This antibody will recognize PDE3B variants. |
| Immunogen Sequence | |



Western blot with Anti-PDE3B Antibody (STJ502269) on PDE3B peptide. Antibody at 1:500 dilution in buffer.