

Anti-PDE3B antibody (1060-1110) {FITC} (STJ502271)

STJ502271

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE3B (1060-1110) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Chicken, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.70 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29516
Gene Symbol	Pde3b
Uniprot ID	PDE3B_RAT
Immunogen	Synthetic peptide taken within amino acid region 1060-1110 on human PDE3B protein.
Immunogen Region	1060-1110
Specificity	This antibody does not cross react with PDE4A variants nor with any other PDE family members. This antibody will recognize PDE3B variants.
Immunogen Sequence	