

Anti-PDE4 antibody {FITC} (STJ502278)

STJ502278

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE4 is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Chicken, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.62 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1,000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide common to most PDE4 subtypes: PDE4A, PDE4D, PDE4B
Immunogen Region	
Specificity	
Immunogen Sequence	