

## Anti-PDE4A antibody (N-Term) {FITC} (STJ502292)

STJ502292

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-PDE4A (N-Term) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.00-1.25 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5141</a>
<b>Gene Symbol</b>	<a href="#">PDE4A</a>
<b>Uniprot ID</b>	<a href="#">PDE4A_HUMAN</a>
<b>Immunogen</b>	N-terminus PDE4A4 synthetic cyclic peptide
<b>Immunogen Region</b>	N-Term
<b>Specificity</b>	Does not cross react with PDE4A1, PDE4A5 or PDE4A8 proteins. Does not cross react with other PDE family members such as PDE4B, PDE4C and PDE4D.
<b>Immunogen Sequence</b>	