

Anti-ENSG00000285188 antibody (C-Term) {FITC} (STJ502319)

STJ502319

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ENSG00000285188 (C-Term) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:2, 500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5143
Gene Symbol	PDE4C
Uniprot ID	PDE4C_HUMAN
Immunogen	Synthetic cyclic carboxy-terminal peptide common to all PDE4C variants.
Immunogen Region	C-Term
Specificity	This antibody detects 3 isoforms of PDE4C subclasses: PDE4C1, PDE4C2, and PDE4C3.
Immunogen Sequence	