

Anti-PDE4D antibody (N-Term) {FITC} (STJ502339)

STJ502339

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE4D (N-Term) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5144
Gene Symbol	PDE4D
Uniprot ID	PDE4D_HUMAN
Immunogen	A 20 amino acid synthetic peptide taken from the N-terminus corresponding to the PDE4D5 variant. A This antibody does not cross-react with other members of the PDE4D family or with other PDE4 proteins (PDE4A, PDE4B and PDE4C) or with other PDE family
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	