

## Anti-PDE6B antibody {Biotin} (STJ502361)

STJ502361

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-PDE6B is suitable for use in ELISA, Confocal Microscopy, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	ELISA, CM, ICC, IF, IHC, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Bovine, Human, Mouse, Pig, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5158</a>
<b>Gene Symbol</b>	<a href="#">PDE6B</a>
<b>Uniprot ID</b>	<a href="#">PDE6B_HUMAN</a>
<b>Immunogen</b>	Synthetic cyclic peptide common to all PDE6 Beta variants
<b>Immunogen Region</b>	
<b>Immunogen Specificity</b>	
<b>Immunogen Sequence</b>	