

Anti-PDE7A antibody {Biotin} (STJ502367)

STJ502367

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE7A is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250 IHC: 1:100-1:250 ICC: 1:100-1:250 IF: 1:100-1:250 CFM: 1:100-1:250
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5150
Gene Symbol	PDE7A
Uniprot ID	PDE7A_HUMAN
Immunogen	Synthetic cyclic peptide common to all PDE7A variants.
Immunogen Region	
Specificity	This antibody is specific for PDE7A variants and will not cross-react with PDE7B or other PDE family members.
Immunogen Sequence	