

## Anti-PDE7A antibody {Biotin} (STJ502369)

STJ502369

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Goat polyclonal antibody anti-PDE7A is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Goat
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.50-0.75 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">18583</a>
<b>Gene Symbol</b>	<a href="#">Pde7a</a>
<b>Uniprot ID</b>	<a href="#">PDE7A_MOUSE</a>
<b>Immunogen</b>	Synthetic cyclic peptide common to all PDE7A variants
<b>Immunogen Region</b>	
<b>Specificity</b>	This antibody is specific for PDE7A variants and will not cross-react with PDE7B or other PDE family members.
<b>Immunogen Sequence</b>	