

Anti-PDE6C antibody (1-75) {Biotin} (STJ502396)

STJ502396

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE6C (1-75) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5146
Gene Symbol	PDE6C
Uniprot ID	PDE6C_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-75 on cone-specific PDE6C protein.
Immunogen Region	1-75
Specificity	
Immunogen Sequence	