

Anti-Phospho-Eph receptor A4 antibody (570-630) {FITC} (STJ502416)

STJ502416

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Eph receptor A4 (570-630) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13838
Gene Symbol	Epha4
Uniprot ID	EPHA4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 570-630 on human EPHA R4 protein.
Immunogen Region	570-630
Specificity	
Immunogen Sequence	