

Anti-Phospho-FGFR2 antibody (520-570) {Biotin} (STJ502427)

STJ502427

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-FGFR2 (520-570) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Cat, Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2263
Gene Symbol	FGFR2
Uniprot ID	FGFR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 520-570 on human Fibroblast Growth Factor receptor 2 protein.
Immunogen	520-570
Region	
Specificity	
Immunogen	
Sequence	