

Anti-Phospho-FGFR4 antibody (530-580) (STJ502432)

STJ502432

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-FGFR4 (530-580) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2264
Gene Symbol	FGFR4
Uniprot ID	FGFR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 530-580 on human Fibroblast growth factor receptor 4 variant.
Immunogen	530-580
Region	
Specificity	
Immunogen	
Sequence	