

Anti-PGAP3 antibody (10-60) (STJ502446)

STJ502446

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PGAP3 (10-60) is suitable for use in Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	DB, ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1,000 DB: 1:50,000 ELISA: 1:50,000 IP: 1:150 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	93210
Gene Symbol	PGAP3
Uniprot ID	PGAP3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 10-60 on human PGAP3 protein.
Immunogen Region	10-60
Specificity	
Immunogen Sequence	



Immunohistochemistry of baboon Colon with Anti-PGAP3 primary antibody (STJ502446) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.