

Anti-Phospho-JAK3 antibody (750-800) (STJ502513)

STJ502513

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-JAK3 (750-800) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3718
Gene Symbol	JAK3
Uniprot ID	JAK3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 750-800 on human JAK3 protein. Phosphorylated tyrosine 785 residue.
Immunogen Region	750-800
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse lung with Anti-Phospho-JAK3 primary antibody (STJ502513) at 1:50 dilution in buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.