

Anti-Phospho-Phospholipase A2 IIA antibody (50-150) {Biotin} (STJ502538)

STJ502538

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Phospholipase A2 IIA (50-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5320
Gene Symbol	PLA2G2A
Uniprot ID	PA2GA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-150 on human Phospholipase A2 membrane associated protein.
Immunogen Region	50-150
Specificity	
Immunogen Sequence	