

Anti-Phospho-MERTK antibody (746-757) (STJ502558)

STJ502558

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-MERTK (746-757) is suitable for use in Confocal Microscopy, ELISA, Flow Cytometry, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, FC, ICC, IF, IHC, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.63 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:750 DB: 1:10, 000 ELISA: 1:10, 000 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65037
Gene Symbol	Mertk
Uniprot ID	MERTK_RAT
Immunogen	Synthetic peptide corresponding to amino acids 746-757 of Human MERTK.
Immunogen Region	746-757
Specificity	
Immunogen Sequence	



Western blot with Anti-Phospho-MERTK Antibody (STJ502558) on mouse macrophage. Antibody at 1:800 dilution.