

Anti-POP4 antibody (140-190) (STJ502585)

STJ502585

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-POP4 (140-190) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Cat, Human, Monkey, Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10775
Gene Symbol	POP4
Uniprot ID	RPP29_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 140-190 on human Ribonuclease P protein subunit p29 protein.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	