

## Anti-Phospho-PDE11A antibody (75-125) {Biotin} (STJ502597)

STJ502597

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-PDE11A (75-125) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.69 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">50940</a>
<b>Gene Symbol</b>	<a href="#">PDE11A</a>
<b>Uniprot ID</b>	<a href="#">PDE11_HUMAN</a>
<b>Immunogen Region</b>	Synthetic peptide taken within amino acid region 75-125 on PDE11A protein. Phosphorylated Ser 117.
<b>Specificity</b>	75-125
<b>Immunogen Sequence</b>	