

Anti-Phospho-PDE5A antibody (120-170) {Biotin} (STJ502622)

STJ502622

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PDE5A (120-170) is suitable for use in ELISA, Immunocytochemistry, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, ICC, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Mouse, Pig, Rat, Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100 CFM: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Phosphorylated synthetic peptide common to all PDE5A variants taken within amino acid region 120-170 on human PDE5A protein.
Immunogen Region	120-170
Specificity	This antibody will detect only PDE5A protein that has been phosphorylated at the PKG site. Will not detect unphosphorylated PDE5A, or other PDE family members.
Immunogen Sequence	