

Anti-Phospho-PGCA antibody (500-550) (STJ502627)

STJ502627

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PGCA (500-550) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4881
Gene Symbol	NPR1
Uniprot ID	ANPRA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human Atrial natriuretic peptide receptor 1 protein.
Immunogen	500-550
Region	
Specificity	
Immunogen	
Sequence	