

Anti-Phospho-RBM4 antibody (300-350) {FITC} (STJ502644)

STJ502644

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-RBM4 (300-350) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5936
Gene Symbol	RBM4
Uniprot ID	RBM4_HUMAN
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 300-350 on human RNA binding motif protein 4.
Immunogen	300-350
Region	
Specificity	
Immunogen	
Sequence	