

Anti-PREX1 antibody (815-836) {FITC} (STJ502650)

STJ502650

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PREX1 (815-836) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57580
Gene Symbol	PREX1
Uniprot ID	PREX1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid region 815-836 of Human PREX1. Peptide was covalently modified.
Immunogen Region	815-836
Specificity	
Immunogen Sequence	