

Anti-PREX2 antibody (720-900) {Biotin} (STJ502652)

STJ502652

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PREX2 (720-900) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.65 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80243
Gene Symbol	PREX2
Uniprot ID	PREX2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 720-900 on PREX2 protein.
Immunogen Region	720-900
Specificity	
Immunogen Sequence	