

Anti-Prominin 1 antibody (700-730) {Biotin} (STJ502665)

STJ502665

GENERAL INFORMATION

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Prominin 1 (700-730) is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA, WB |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|-----------------------------------------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.50-0.95 µg/µl |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500-1:1, 000 DB: 1:2, 000 ELISA: 1:2, 000 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------------------------------------------------------|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Synthetic peptide taken within amino acid regions 700-730 on CD133 protein |
| Immunogen Region | 700-730 |
| Specificity | |
| Immunogen Sequence | |