

Anti-Brk antibody (400-450) (STJ502711)

STJ502711

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Brk (400-450) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	20459
Gene Symbol	Ptk6
Uniprot ID	PTK6_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 400-450 on mouse Protein-Tyrosine Kinase 6 isoform.
Immunogen	400-450
Region	
Specificity	
Immunogen	
Sequence	