

Anti-Phospho-TICAM2 antibody {Biotin} (STJ502724)

STJ502724

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-TICAM2 is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:750 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100302736/353376
Gene Symbol	TICAM2
Uniprot ID	TCAM2_HUMAN
Immunogen	Phosphorylated synthetic peptide corresponding to unique amino acid sequences on TRIF-related adapter molecule protein.
Immunogen Region	
Specificity	
Immunogen Sequence	