

Anti-QSOX1 antibody (500-550) (STJ502741)

STJ502741

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-QSOX1 (500-550) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5768
Gene Symbol	QSOX1
Uniprot ID	QSOX1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-550 on human QSOX1 protein.
Immunogen Region	500-550
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat small intestine with Anti-QSOX1 primary antibody (STJ502741) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.