

## Anti-Dexas1 antibody (256-278) (STJ502750)

STJ502750

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Dexas1 (256-278) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Bovine, Human, Monkey, Panda

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.55 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	19416
<b>Gene Symbol</b>	Rasd1
<b>Uniprot ID</b>	RASD1_MOUSE
<b>Immunogen</b>	Synthetic peptide: HSDLMYIREKTSVGSQAKDKERC corresponding to amino acids 256-278 of Mouse Dexas 1 protein.
<b>Immunogen Region</b>	256-278
<b>Specificity</b>	Does not cross react to its closely resembling RHES protein.
<b>Immunogen Sequence</b>	



Western blot with Anti-Dexas1 Antibody (STJ502750) with RASD1 peptide. Antibody at 1:1000 dilution in buffer.