

Anti-Dexas1 antibody (256-278) {Biotin} (STJ502751)

STJ502751

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dexas1 (256-278) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey, Panda

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	19416
Gene Symbol	Rasd1
Uniprot ID	RASD1_MOUSE
Immunogen	Synthetic peptide: HSDLMYIREKTSVGSQAKDKERC corresponding to amino acids 256-278 of mouse Dexas 1 protein.
Immunogen Region	256-278
Specificity	Does not cross react to its closely resembling RHES protein.
Immunogen Sequence	