

Anti-RBP2 antibody (200-250) {FITC} (STJ502758)

STJ502758

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RBP2 (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74152
Gene Symbol	Stra6l
Uniprot ID	STR6L_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 200-250 on mouse Retinol-Binding Protein 4 Receptor 2/RBP4 receptor 2/STRA6-Like Protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	