

Anti-REST antibody (1000-1097) (STJ502771)

STJ502771

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-REST (1000-1097) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-0.54 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5978
Gene Symbol	REST
Uniprot ID	REST_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1000-1097 on human REST protein.
Immunogen Region	1000-1097
Specificity	
Immunogen Sequence	



Western blot with Anti-REST Antibody (STJ502771) with REST peptide. Antibody at 1:500 dilution in buffer.