

Anti-RIC8A antibody (500-535) {FITC} (STJ502779)

STJ502779

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIC8A (500-535) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.60-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:250 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60626
Gene Symbol	RIC8A
Uniprot ID	RIC8A_HUMAN
Immunogen	Synthetic peptide taken within amino acids 500-535 on RIC8A protein
Immunogen Region	500-535
Specificity	
Immunogen Sequence	