

Anti-RICH2 antibody (1-60) (STJ502780)

STJ502780

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RICH2 (1-60) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:250-1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	216831
Gene Symbol	Arhgap44
Uniprot ID	RHG44_MOUSE
Immunogen	Synthetic peptide corresponding to N terminal amino acids 1-60 of Mouse RICH2.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	



Western blot with Anti-RICH2 Antibody (STJ502780) with 1) Recombinant RICH1, 2) Recombinant RICH2, 3) CaCO2 cells, 4) Jeg3 cells. Antibody at 1:500 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081