

Anti-RICH2 antibody (1-60) {Biotin} (STJ502781)

STJ502781

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RICH2 (1-60) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:250-1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	216831
Gene Symbol	Arhgap44
Uniprot ID	RHG44_MOUSE
Immunogen	Synthetic peptide corresponding to N terminal amino acids 1-60 of Mouse RICH2.
Immunogen Region	1-60
Specificity	
Immunogen Sequence	