

Anti-Robo3 antibody {Biotin} (STJ502784)

STJ502784

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Robo3 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Monkey, Panda

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.52-0.63 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29546
Gene Symbol	Homer1
Uniprot ID	HOME1_RAT
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human RIG1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	