

Anti-RIOK1 antibody (450-500) (STJ502786)

STJ502786

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RIOK1 (450-500) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83732
Gene Symbol	RIOK1
Uniprot ID	RIOK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-500 on human Serine/threonine-protein kinase RIO1 protein.
Immunogen Region	450-500
Specificity	
Immunogen Sequence	