

## Anti-RIP140 antibody (50-100) {Biotin} (STJ502794)

STJ502794

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-RIP140 (50-100) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IHC, IF, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.60 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">268903</a>
<b>Gene Symbol</b>	<a href="#">Nrip1</a>
<b>Uniprot ID</b>	<a href="#">NRIP1_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 50-100, 740-690 or 1050-2000 on mouse receptor interacting protein 140.
<b>Immunogen Region</b>	50-100
<b>Specificity</b>	
<b>Immunogen Sequence</b>	